




The Balfour retractor is a self-retaining abdominal retractor commonly used in abdominal and pelvic surgeries. It consists of blades that are inserted into the abdominal incision and then spread apart to hold back the abdominal wall, providing a clear view and access to the internal organs. The Balfour retractor helps maintain exposure during procedures like laparotomies, allowing surgeons to work effectively within the abdominal cavity.



Blades

Product Number	Dimensions
TS/BFR/BL/01	40 x 80 mm
TS/BFR/BL/02	25 x 80 mm
TS/BFR/BL/03	55 x 80 mm
TS/BFR/BL/04	60 x 65 mm

*Blade size can be custom made as per requirement.